



# COAX

SEMI-AUTOMATED COAGULOMETERS

# SEMI-AUTOMATED COAGULOMETERS



## With 1,2 or 4 optical channel

- Prepared for the daily routine and upcoming requirements.
- High quality in the results.
- Nearly maintenance free

## COAX

Coagulation is a change of physical state of the blood due to the conversion of a soluble plasma protein, fibrinogen, into a solid gel, fibrin.

The management and control of anticoagulant therapy and the assessment of pre and post surgical states, among others requires a proper evaluation of the coagulation cascade.

Several tests help the physician in the diagnosis of alterate coagulation states and management of coagulopathy.

The coagulation reagents have been specifically validated to BioSystems coagulometers.

## Prothrombin Time (PT)

|    | Presentation | Code  |
|----|--------------|-------|
| PT | 4x5 mL       | 61001 |




### Principle of the method:

The addition of calcium thromboplastin to plasma induces the formation of the fibrin clot. The method measures the clot formation time.

### Intended use:

- Screening assay used to monitor oral anticoagulant therapy.
- It helps detect and diagnose a bleeding disorder.

### Procedure




## Activated Partial Thromboplastin Time (APTT)

|                             | Presentation             | Code  |
|-----------------------------|--------------------------|-------|
| APTT A                      | 4x4 mL                   | 61004 |
| APTT B (CaCl <sub>2</sub> ) | 4x16 mL                  | 61005 |
| APTT                        | A (4x4 mL) + B (1x16 mL) | 61009 |



### Procedure



### Principle of the method:

The addition of the phospholipid cephalin to plasma samples in the presence of calcium and an activator induces the formation of the fibrin clot. The method measures the clot formation time.

### Intended use:

- Screening assay used in the monitoring of heparin therapy.
- As part of investigation of a possible bleeding disorder.

## Activated Partial Thromboplastin Time (APTT)

|                  | Presentation             | Code  |
|------------------|--------------------------|-------|
| FIB A            | 4x2 mL                   | 61002 |
| FIB B (Imidazol) | 4x15 mL                  | 61003 |
| FIB              | A (4x2 mL) + B (4x15 mL) | 61004 |



### Procedure



### Principle of the method:

The addition of calcium thromboplastin to plasma induces the formation of the fibrin clot. The method measures the clot formation time.

### Intended use:


- Screening assay used to monitor oral anticoagulant therapy.
- It helps detect and diagnose a bleeding disorder.

## Thrombin Time (TT)

|    | Presentation | Code  |
|----|--------------|-------|
| TT | 4x3 mL       | 61000 |



### Procedure



### Principle of the method:

Addition of human thrombin to plasma samples induces the formation of fibrin clot. The method measures the clot formation time.

### Intended use:

- To evaluate the level and function of fibrinogen.
- To detect heparin contamination.
- As part of investigation of a bleeding or thrombotic episode.

## Calibrators and Controls

|                  | Presentation | Code  |
|------------------|--------------|-------|
| Calibrator       | 4x1 mL       | 61006 |
| Calibrator       | 2x1 mL       | 61015 |
| Control I        | 4x1 mL       | 61007 |
| Control II       | 4x1 mL       | 61008 |
| Control level I  | 2x1 mL       | 61012 |
| Control level II | 2x1 mL       | 61013 |

The Coagulation Calibrator is a lyophilized pooled human plasma containing component concentrations suitable for the calibration of measurement procedures.

The Coagulation Control is a lyophilized human plasma with stabilizer suitable for the quality control of the clinical laboratories. The product is intended for intralaboratory quality control purposes only and is supplied with intervals of suggested acceptable values.



## Specifications

| Code                    | 85001                      | 85002              | 85004             |
|-------------------------|----------------------------|--------------------|-------------------|
| Optical channels        | 1                          | 2                  | 4                 |
| Wavelength (nm)         | 620 (red)                  | 405 (UV)           | 405 (UV)          |
| Global Coag. Tests      | PT, APTT, TT, FIB          | PT, APTT, TT, FIB  | PT, APTT, TT, FIB |
| Specific Coag. Tests    | -                          | individual factors |                   |
| Chromogenic Coag. Tests | -                          | AT,PC              |                   |
| Latex based tests       | D-Dimer                    |                    |                   |
| Display                 | Color touch screen display |                    |                   |
| Dimensions              | 230 x 140 x 90 mm (l,b,h)  |                    |                   |
| Interfaces: RS 232 (2x) | Printer, Barcode reader    |                    |                   |
| USB (2x)                | Network, Firmware update   |                    |                   |

## Consumables

| Product             | Code  |
|---------------------|-------|
| 1 pack 500 cuvettes | 85020 |

### MOGADISHU OFFICE

+252 61 5000491  
+252 1 656503

### HARGEISA OFFICE

+252 63 4078539  
+252 2 512422

### GAROWE OFFICE

+252 90 7742425  
+252 90 7662918

### DJIBOUTI OFFICE

+253 77 021504  
+253 77 189595