

# **Easy Stat**

Blood Gas Analyzer



# Easy inside and out

Medica's EasyStat analyzer measurespH,
PCO2, PO2, Na+, K+, Ca++or Cl-andHct,
and calculates additional parame-ters.
Patient parameters, including
FIO2, patient ID, patient temperature,
%FIO2, drawing site and other information
can be entered using the digital keypad and
integrated with patient results.
Measuredand calculated results are displayed andprinted•]

EasyStat focuses on the laboratory's need to deliver sample results economi-cally and efficiently•

The sophistication and performancerequired by today's busy, demanding health care environment have been packaged in a new compact formatwith a small footprint to save space•

Liquid calibrants are packaged in a convenient reagent module, eliminatinggas tanks• All components are combined intothree simple modules, easily accessi-ble by the user.

Routine maintenance islimited to the replacement of electrodes and a single pump tube.•

Simple menus guide the user through analyzer operation.

Unique electrode design with no membranes to change, combined with areagent module with over 1,000-sample capacity ensure economical operationand low cost per sample



HOME MENU	DIAGNOSTICS	SETUP MENU
	2	
1. ANALYZE SAMLE	1. TEST COMPONENTS	1. USER OPTIONS
2. ANALYZE QC	2. TEST FLUIDICS	2. CONFIGURATION
3. CALIBRATE	3. SENSOR STATUS	3. PRINTER OPTIONS
4. DAILY CLEANER	4. PRIME FLUIDS	4. SET REFERENCE LIMITS
5. SECOND MENU	5. PRINT mV's	5. PATIENT INFORMATION
		6. DELETE DATA






# Flexible sofware

flexible easy stat software can be programed to coform with the lab's established practices. The sofware allows selection of desired option,including:

Language Reference limits Correlation with other methods Parameters Reported

Quanlity control limits operation ID Patient data



Electrode design permits fast,fail-safe installation

# EasyStat can be used and maintained by anyone, anytime, anywhere•

selects calibrants and rinse solution

Removal of the three plug-in modules -Reagent Module, Sensor Module and ValveModule —is accomplished without tools.

Medica's integral membrane design means that electrodes are disposable and require no maintenance. Electrodes snap in and out for easy replacement.

Sample path has been simplified—only one length of pump tubing requires periodic replacement..

Innovative design simplifies maintenance, addressing the needs of the remote laboratory which has limited access to technical service personnel.

Diagnostic software displays component status, assuring quick troubleshooting. Modularitymakes assembly and disassembly quick and easy. There is no need for expensive service contracts.



## Comprehensive management of patient, quality control and maintenance data.

The EasyStat quality control program calculates and stores complete statistics for thelast 31 days of quality control results at each of three levels.

A printed Levey-Jenningschart visually identifies trends.•

The data management program compares all patient results with ranges stored in memoryand flags out-of-range results.

Results are stored in memory for up to 64 patients.•

Exchange of components, calibration and other events are documented for regulatorycompliance.•

Bar code scanner option permits rapid, accurate input of patient, operator and QC data.

# Easy outside

### Blood gas analyzer operation has never been simpler

The Universal Sampler adapts to both syringe and capillary samples without adaptors.

The sample probe's self-wiping feature provides convenience, sample integrity and user safety.

The simple, yes/no user-prompting menu makes rapid training of new personnel possible.

# Compact Reagent Module for convenience, economy and safety

Bulky gas tanks are replaced with liquid, tonometered calibrants, packaged in a convenientReagent Module that also collects waste, protecting the user from biological hazards. The EasyStatautomatically tracks date code and calibrant usage with the Reagent Module's solid state memory. Operation without interruption is assured.

#### **Portable**

Light weight design with optional handle allows use in any setting. Blood gas analyzer operation has never been simpler











# Specifications

CLIA Classification: Moderate complexity

Sample Type: Whole blood

Sample Size: 120 µL Syringe mode/95 µL Capillary mode

### **Measured Parameters and Limits**

PO2 5- 700 mmHg
PCO2 5.0- 150.0 mmHg
pH 6.500- 8.000 pH units
Hct 10- 70%

Hct 10-70% Na+ 80-200 mmol/L K+ 1.0-20.0 mmol/L

Ca++ 0.25- 5.00 mmol/L CI- 50.0- 150.0 mmol/L

### **Calculated Parameters**

THb (Total Hemoglobin) 3.3-23.3 g/d

LpH (T) (pH temperature corrected)

PCO2(T) (PCO2temperature corrected)

PO2(T) (PO2temperature corrected)

TCO2 (Total Carbon dioxide) 0– 50 mmol/L HCO3- (Bicarbonate) 0– 50 mmol/L

BEb (Base Excess in blood) -25.0 to 25.0 mmol/L
BEecf(Base Excess in extracellular fluid) -25.0 to 25.0 mmol/L

SBC (Standard Bicarbonate) 0- 50 mmol/L

%SO2c(Oxygen Saturation) 40.0-100.0% (calculated at normalP50)

CtO2 (Oxygen Content) 3.0-30.0 mL/dL A-aDO2 (Alveolar arterial oxygen gradient) 0-700 mmHg

RI (Respiratory Index) 0.0-70.0

Ca++(7.4) (for 7.2<pH<7.6) 0.22-5.58 mmol/L

#### **Input Parameters**

Patient Temperature (20- 45°C)
Hemoglobin (3.0- 30.0 g/dL)
FIO2 (Fraction Inspired Oxygen) (10- 100%)
Patient ID (14 digits)
Operator ID (14 digits)
Time Drawn (00:00)

Sample Source (arterial, mixed venous, venous)
Sample Type (radial, brachial, femoral, arterial line)

### **BLOOD GAS ANALYZER**



Sample Temperature Control: 37.0°C ± 0.2°C

Ambient Conditions: 15–30°C (59–86°F), 500–800 mmHg (max 15 PSI)

5–85% relative humidity, non-condensing atmospheric

air environ ment (21% O2)

Analysis Time: <120 seconds

Data Storage: 64 Patient results with Operator ID, Patient ID,

Date and Time QC—up to 93 results for each Level

(Blood Gas/Electrolytes 1, 2, 3, Hct 1, 2)

Calibration: Automatic or On-Demand

Input/Output: Numeric keypad, graphic display, 27 column thermal line

printer, barcode reader port, RS-232 computer interface

port

Power: 100/115~VAC, 50-60 Hz, 0.8 A or 220~VAC, 50-60 Hz, 0.4

A Refer to the chassis serial number label

for the voltage tha has been factory set on your analyzer, and for proper fuse re place ment.

Size & Weight: 14.5" W x 12.5" H x 7.0" D (37cm W x 32cm H x 18cm D),

17lbs (7.7 kg) with Reagent Module

